

# OASYS with Invisalign Treatment



- The Dentist / Orthodontist will take Impressions with the Invisalign Trays in Place
- The Upper Invisalign Tray will take the place of the Upper Oasys Splint
- The lab will blockout and fabricate another Invisible Retainer over the Lower Model
- The OASYS will be fabricated. Retentive Clasps will be included for some Retention on the Lower.
- The Dentist / Orthodontist will use PVS / Bite material to Add Retention.
- If ThermoFlex Splint Material is used, Microwave Water can soften splint, and the splint squeezed to add some retention
- The Trays will be changed per Treatment Protocol
- The Lower Buccal & Lingual of the Splint can be trimmed, as needed.
- More PVS / Bite material is added for retention.
- The Lower Splint will be Relined to House the Finished Invisalign Tray or the Lower OASYS Splint will be Relined to the New Dentition Alignment for Sleep Apnea Treatment and Lower Arch Retention.