



FDA 510K #971794

The **EMA Devices** fabricated by **Dream Systems** use the **Myerson Hooks and Straps**. The standard appliance design has been modified, incorporating a flat plane design of the splints for added strength and function, The **ThermoFlex Hard Material** is used upon request and in circumstances where there are impression flaws, edentulous areas and where claspings can be of benefit. **Hard Acrylic with Clasps** is also available for Digital Printed Models. An **Intro Strap Pack** is included with each case, which will provide **5mm of additional advancement** from the **Protrusive Bite** provided. The Strap Sequence and Color Flexibility of the Straps will vary by individual case. Straps may be included to decrease the protrusion or strain of a more rigid strap. Additional Straps can be ordered as determined by effectiveness.

Strap Sizes

- 21mm
- 20mm
- 19mm
- 18mm
- 17mm
- 16mm
- 15mm
- 14mm
- 13mm

White—Very Flexible—Provided Only Upon Request

Yellow—Less Flexible—Provided if the Protrusive Bite taken is advance of end to end and also used to decrease in protrusive strain.

Blue—Standard Strap Flexible—Included for 5mm extension.

Clear—Rigid— Provided Only Upon Request for Less Movement Side to Side and Longevity of the Straps

STRAPS



Vacuform Splints

- Standard
- 2.0 MM
- Bruxer
- 2.5 MM
- HardSoft
- 2.0 MM
- Digital or or
- Good Detail in the Impressions
- Required



EMA First Step 90 Day Trial Device

Includes Three pairs of Yellow Straps for advancement for an additional 4mm.



* A Comfortable Protrusive Bite Record with a Minimum 4mm Incisal Opening is recommended.

Hard Acrylic Splints with Clasps

is available for Digital Printed Casts

ThermoFlex Hard Splints

Ball or Arrow Clasps Included, More Ideal for Missing Teeth and More Severe Undercuts,

660 Commerce Drive, Ste. C
 Roseville, CA 95678
 (916) 865-4528 (916) 872-9569
www.dreamsystemsdentallab.com