


The **Dream Systems KAVA Dorsal™ Appliance** is **FDA 510K Cleared # K182661**, It incorporates Metal Reinforce Dorsal Fins and an Angled Advancement Screw, with Advancement to **7mm**. They hold the lower jaw forward during sleep, but allows the patient to open and close his/her mouth and have some side to side movement This will help to open the airway in back of the throat, while pulling the tongue forward.

- METAL REINFORCED FINS
- 7MM ADVANCEMENT SCREW
- RIGHT & LEFT SIDE ACTIVATION UPWARDS
- 2MM SIDE TO SIDE MOVEMENT



- 7MM ADVANCEMENT SCREW
- 10 KEY HOLE TURNS FOR 1MM ADVANCEMENT



NASAL DILATORS CAN BE ADDED FOR IMPROVED NASAL BREATHING

